


Bacteriophage Rapid Anthrax Viability Evaluation System

BRAVES


Patented LumiPhage™ Technology



Light it up!

- Minimal training & hands on time (<10 min)
- High-throughput (96 samples ≤2 h)
- Detects 10¹ CFU in ≤8 h
- Functions with 'dirty' samples
- Detects viable *B. anthracis* pre- and post-decontamination
- Tested with steel, glass, plastic, cement, wood, water & soil

Simple workflow:



Rapid results: 10¹ CFU in ≤8 h

